



## NITRATE/NITRITE KIT

OCTA-SLIDE 2, 0.25-10.0 ppm

CODE 3519-01


QUANTITY	CONTENTS	CODE
2 x 120 mL	*Mixed Acid Reagent	*V-6278-J
5 g	*Nitrate Reducing Reagent	*V-6279-C
5 g	*Color Developing Reagent	*V-6281-C
2	Spoons, 0.1g, plastic	0699
3	Test Tubes, 2.5-5-10 mL, w/caps, plastic	0106
1	Water Sample Bottle	0688
1	Dispenser Cap	0692
1	Nitrate-N Octa-Slide 2 Bar, 0.25-10.0 ppm	3109-01
1	Octa-Slide 2 Viewer	1101

**\*WARNING:** Reagents marked with an \* are considered to be potential health hazards. To view or print a Safety Data Sheet (SDS) for these reagents go to [www.lamotte.com](http://www.lamotte.com). Search for the four digit reagent code number listed on the reagent label, in the contents list or in the test procedures. Omit any letter that follows or precedes the four digit code number. For example, if the code is 4450WT-H, search 4450. To obtain a printed copy, contact LaMotte by email, phone or fax.

Emergency information for all LaMotte reagents is available from Chem-Tel: (US, 1-800-255-3924) (International, call collect, 813-248-0585)

To order individual reagents or test kit components, use the specified code number.

### USE OF THE OCTA-SLIDE 2 VIEWER



The Octa-Slide 2 Viewer should be held so non-direct light enters through the back of the Viewer. Slide the Octa-Slide 2 Bar into the Viewer. Insert the reacted sample into the top of the Viewer. Match the color of the reaction to the color standards.

**WARNING!** This set contains chemicals that may be harmful if misused. Read cautions on individual containers carefully. Not to be used by children except under adult supervision.

Nitrites can cause serious interference in the Nitrate test. To obtain the actual nitrate-nitrogen concentration, follow Procedures A & B. Subtract the nitrite reading obtained in Procedure B (as ppm nitrate-nitrogen) from the nitrate reading obtained in Procedure A (as nitrate-nitrogen).

The best results are obtained when the sample and reagents are at  $23 \pm 2^{\circ}\text{C}$ .

### **PROCEDURE A - NITRATE TEST**

---

1. Insert the Nitrate-N Octa-Slide 2 Bar (3109-01) into the Octa-Slide 2 Viewer (1101).
2. Fill sample bottle (0688) with sample water.
3. Fill test tube (0106) to the 2.5 mL line with water from the sample bottle.
4. Dilute to the 5.0 mL line with \*Mixed Acid Reagent (V-6278). Cap and mix. Wait two minutes.
5. Use the 0.1 g spoon (0699) to add one level measure (avoid any excess) of \*Nitrate Reducing Reagent (V-6279). Cap and invert the tube 50-60 times in one minute. Wait 10 minutes.
6. Insert the test tube into the Octa-Slide 2 Viewer.
7. Match sample color to a color standard. Record as ppm Nitrate-Nitrogen.
8. To convert to ppm Nitrate ( $\text{NO}_3$ ), multiply the test result by 4.4.

NOTE: To obtain the actual nitrate concentration, follow with Procedure B and subtract the nitrite reading from the reading obtained in the nitrate test.

NOTE: Thoroughly clean and rinse test tubes after each use.

### **PROCEDURE B - NITRITE TEST**

---

1. Insert the Nitrate-N Octa-Slide 2 Bar (3109-01) into the Octa-Slide 2 Viewer (1101).
2. Fill the sample bottle (0688) with sample water.
3. Fill test tube (0106) the 2.5 mL line with water from the sample bottle.
4. Dilute to the 5.0 mL line with \*Mixed Acid Reagent (V-6278). Cap and mix.
5. Use the 0.1 g spoon (0699) to add one level measure (avoid any excess) of \*Color Developing Reagent (V-6281). Cap and invert the tube 50-60 times in one minute. Wait 10 minutes.
6. Insert the test tube into the Octa-Slide 2 Viewer.
7. Match sample color to a color standard. Record as ppm Nitrite as Nitrate-Nitrogen.
8. Use the following table to convert reading to Nitrite-Nitrogen. Record as ppm Nitrite-Nitrogen.

<b>Comparator Reading</b>	<b>Divide By</b>
<b>0.5 or less</b>	<b>3.5</b>
<b>0.5-1.5</b>	<b>4.6</b>
<b>over 2.0</b>	<b>5.8</b>

NOTE: If the nitrate-nitrogen reading is higher than 10.0 the sample can be diluted 1:1 with nitrate and nitrite free water, the test run and the resultant value obtained multiplied by 2.

NOTE: Thoroughly clean and rinse test tubes after each use.

### **LaMOTTE COMPANY**

Helping People Solve Analytical Challenges®

PO Box 329 · Chestertown · Maryland · 21620 · USA  
800-344-3100 · 410-778-3100 [Outside U.S.A.] · Fax 410-778-6394

Visit us on the web at [www.lamotte.com](http://www.lamotte.com)