

Code 3038-G

25 Test Strips

5-WAY
Nitrate, Nitrite, pH, Alkalinity and Total Hardness



Insta-Test[®] ANALYTIC

Fresh & Salt	0	20	40	80	160	200	
	OK	OK	OK	OK	OK	OK	
Fresh & Salt	0	0.5	1	3	5	10	
	OK	OK	OK	OK	OK	OK	
Fresh	6.0	6.5	7.0	7.5	8.0	8.5	9.0
	OK	OK	OK	OK	OK	OK	OK
Salt	6.0	6.5	7.0	7.5	8.0	8.5	9.0
	OK	OK	OK	OK	OK	OK	OK
Fresh & Salt	0	40	80	120	180	240	
	OK	OK	OK	OK	OK	OK	
Fresh	0	30	60	120	180		
	OK	OK	OK	OK	OK		

NO₃
Nitrate
ppm

NO₂
Nitrite
ppm

pH

Alk
Alkalinity
ppm

TH
Hardness
ppm

- Immerse strip. SWIRL 2 TIMES.
Remove with pads face up.
DO NOT SHAKE OFF EXCESS WATER. Immediately read Hardness > Alkalinity.
Wait 30 seconds, read pH > Nitrite > Nitrate.
- Keep wet fingers out of the vial
 - Close vial tightly after removing strip
 - Store in a cool, dry place

LaMotte Company
Chesterstown, Maryland 21620 USA
www.lamotte.com



Drop holding box and dotted tolerance line